

Application Note

Laser Mode Analysis with Photon Beam Profilers

How Gaussian is Your Beam?

The Gaussian form $I = e^{-2(r/w)^2}$ gives the intensity at each value of r . The parameters are illustrated in Figure 1: r is the radius at a particular intensity I ; w is the radius at the $1/e^2$ intensity.

One can measure the beam diameter at two intensity values and obtain their ratio. The ratio can then be compared to that of a perfect Gaussian beam.

If, for example, we select the $1/e^2$ and 53.5% intensity values, the theoretical beam diameter ratio for a Gaussian beam is found by solving the equation for $I = 0.535$ (the intensity clip level setting):

$$\begin{aligned} \ln(0.535) &= \ln e^{-2(r/w)^2} \\ -0.62548 &= -2(r/w)^2 \\ 0.3127 &= (r/w)^2 \\ 0.5592 &= (r/w) \\ R = \frac{2r}{2w} &= \frac{D(53.5\%)}{D(13.5\%)} = 0.559 \\ \ln(x) &= \log \text{ to base } e \text{ of quantity } x \end{aligned}$$

For profiles close to a pillbox shape, the measured ratio will be greater than 0.559; for waveforms close to a triangle shape, the ratio will be less than 0.559. You can calculate any ratio set you choose by solving the more general equation: $I = e^{-a^2}$ at the two clip levels of interest and determining the ratio r . For illustration, let's assume we'd like to use the 90% and 10% intensity beam diameters; solve for a :

$$\begin{aligned} \ln I &= -a^2 \\ a^2 &= -\ln I \\ a &= \sqrt{-\ln I} \end{aligned}$$

Substituting $I = 90\%$

$$\begin{aligned} a(90\%) &= \sqrt{-\ln(0.90)} \\ a(90\%) &= 0.3245 \end{aligned}$$

With $I = 10\%$:

$$\begin{aligned} a(10\%) &= \sqrt{-\ln(0.10)} \\ a(10\%) &= 1.5174 \end{aligned}$$

The ideal Gaussian ratio is then:

$$R = \frac{a(90\%)}{a(10\%)} = \frac{0.3245}{1.5174} = 0.2138$$

For some applications of lasers, it is very important to know the mode mixture; a Photon beam profiler can be used to detect this mixture. To illustrate, we calculated the ideal profile, and then added increments of TEM_{01} mode to the TEM_{00} beam.

Figures 2-5 illustrate profiles with decreasing TEM_{01} content; the beam diameter decreases and the slopes change; both these changes are detectable with a Photon profiler. Figure 6 illustrates a two-dimensional intensity contour map. Both the pure Gaussian and the Gaussian +5% TEM_{01} contours are shown. Each contour line represents a 5% change in intensity. Visually, it is easy to see that the beam is elliptic; and practically, that the ellipticity can be measured by a Photon profiler.

Repeatability of beam diameter measurements for a standard BeamScan with a 2.5mW commercial HeNe laser is $\pm 3\mu\text{m}$ in $680\mu\text{m}$, or approximately 1/2% at $1/e^2$ diameter. Many parameters can degrade this performance, but let's assume this repeatability and ask what change in Gaussian ratio, R , is significant? Assume the test is made at the 53.5% and $1/e^2$ intensity. The ideal ratio is 0.559; significant ratios are less than 0.551 and greater than 0.567.




Figure 1
Pure TEM_{00} and 99% TEM_{00}
Plus 1% TEM_{01}




Figure 2
Pure TEM00 and 60% TEM₀₀
Plus 40% TEM₀₁




Figure 3
Pure TEM00 and 80% TEM₀₀
Plus 20% TEM₀₁




Figure 4
Pure TEM00 and 90% TEM₀₀
Plus 10% TEM₀₁




Figure 5
Pure TEM00 and 95% TEM₀₀
Plus 5% TEM₀₁

Figure 6
Pure TEM₀₀ and 95% TEM₀₀ plus
5% TEM₀₁ Intensity Contours
(5% intensity change/contour line)

