



## Manual Measuring Stage MTM 150 × 150 MR with Integrated Measuring System

### Specifications

Measuring range	150 × 150 mm (6" × 6")
Accuracy	±1 µm
Resolution	0.05 µm (smallest step size)
Drive	coaxial drive, rack and pinion with patented coarse adjustment
Material	aluminium
Surface	anodic coating, black lacquered
Load capacity	3.0 kg
Weight	4.8 kg

### Order Information

Order No.: 00-31-607-0000

### Product Features

- ▶ highest accuracy through integrated linear measuring systems
- ▶ for reflected and transmitted light applications
- ▶ equipped with patented coarse adjustment
- ▶ incl. digital readout unit PROFILER SCD

### Accessories

Glass plate (included in delivery)  
 Order No.: 00-24-601-0810

Digital readout unit PROFILER SCD  
 (included in delivery)  
 Order No.: 00-76-900-4830

For further accessories please visit:  
[www.marzhauser.com](http://www.marzhauser.com)

### Dimensions

