

Motorized Microscope Stage EK 14 mot for Meiji MT4000–MT8500

Specifications

Travel range	75 × 30 mm
Travel speed	max. 120 mm/s
Repeatability	< 5 µm (unidirectional) ¹
Accuracy	±20 µm ¹
Orthogonality	< 10 arcsec
Motor type	2-phase stepper motor
Drive	rack – pinion
Limit switches	integrated Hall sensors
Material / surface	aluminium / anodic coating, black lacquered
Worktop	wear-free CERASIST® coating
Weight	approx. 2.5 kg

Order Information

Order No.: 96-27-322-0000

Product Features

- ▶ CERASIST® coating, free of wear
- ▶ travel range limitation by integrated Hall sensors
- ▶ optionally integrable coaxial drive for ergonomical handling (adjustable in height and friction)
- ▶ suitable for transmitted and reflected light applications


Accessories

Object holder for 1 slide 76 × 26 mm (included in delivery)

Order No.: 00-27-322-2000

Electrical Coaxial Drive

Order No.: 00-27-322-1000



For more accessories for microscopy please visit: www.marzhauser.com