28 YEARS OF MAKING CRYOGENICS FOR RESEARCHERS www.cryoindustries.com

CRYOCOOL-PC



With transfer line for remote cooling

High Efficiency Nitrogen CRYOCOOL-PC Cold Stream Transfer System:

Small personal size with Never ice 'warm-tip' (no shield gas needed - ever) - ice free nozzle without using shield gas

COLD constant flow gas stream at 95K or less

or carry dewar unit with integral x-y-z stage

Temperature stability better than +/- 0.1K (C)

Large nozzle diameter of 0.354 inch (~9 mm) inside diameter 0.312 inch (~ 8 mm)

Vacuum and superinsulated for highest efficiency

All welded stainless steel transfer line for long term reliability

Charcoal vacuum maintenance cryopump in insert leg

Compact Tabletop Design!



Liquid Nitrogen LED level monitor





Modular design for easy cleaning and service

