


Non-Variable Temperature Optical Pour Fill Liquid Nitrogen Cryostat Sample in Vacuum



System Features:

- 77K base operating temperature
- 0.4 liter liquid nitrogen capacity
- Pour fill the LN₂ reservoir - virtually unlimited runtime
- Compact hand-held style
- Vacuum shroud with 3.25" square bottom section with (4) optical 1.62 inch diameter clear view CFQ suprasil window ports
- Instrumentation header with:
 - (1) 10-pin electrical feedthroughs
 - (3) Blank feedthrough ports
 - Safety pressure relief
 - Evacuation valve with NS25 (ISO-KF 25) interface
 - Copper optical 'T' style holder with 5mm diameter hole
 - 2.56" diameter sample mount with provision for temperature sensor
- Complete system test
- One year warranty